

Article code	Z072910F
Article name acc. DoC	MS 10
Class acc. Directive 93/42/EEC	IIa

Article code	Article name	Material
C11503	Solvent	●THF (Tetrahydrofuran), also 1,4-Epoxy butane
C30910	Tube 2.94 x 4.09 mm	●PVC (Polyvinylchloride), with plasticizer DEHP
C31514	Spring clamp	●Housing: POM (Polyoxymethylene) ●Colourant ●Spring: stainless steel
J1010110		
C10178	Drip chamber upper part 21 E	●PS (polystyrene) ●Colourant white, 50 % TiO ₂ in HIPS (high impact polystyrene)
C10221	moulded ring white for drip chamber 21	●PS (Polystyrene) ●Colourant, carrier material PS (Polystyrene), contains 50% TiO ₂
C10265	Drip chamber lower part 21 V92	●SBS (Styrene Butadiene Copolymer)
C10345	air ventilation flap	●1-propene homopolymer: PP (Polypropylene) ●1-propene copolymer (heterophasic copolymer)
C10539	Protective cap	●PP (Polypropylene)
C10808	Filter 15 µm drip chamber 21	●Filter: PA 6.6 (Polyamide) = Nylon Filter ●housing: PA 6 (Polyamide)
C10836	ventilation filter 3µm	●Membrane Type: Hydrophobic ●Material Type: ePTFE (expanded Polytetrafluorethylene) Membrane ●Membrane Support Type: Non-woven Polyester
C11504	Solvent / Lubricant	●Isopropyle alcohol
C11511	Lubricant	●Aminofunctional Siloxane
J1010810		
C10405	roller clamp housing, orange	●PS (polystyrene) ●PS (Polystyrene) with silicone slip-masterbatch; ●masterbatch orange
C10412	Roller clamp wheel	●PS (polystyrene); ●PS (Polystyrene) with silicone slip-masterbatch
J1010914		
C10564	Protective cap LL pos. with 1.2 µm filter	●Housing: PP (Polypropylene) + Masterbatch ●Filter: Haemophobic/hydrophobic membrane; 1,2 µm PVDF (Polyvinylidene difluoride) Fortex Membrane
C30605	Revolving Luer Lock	●luer cone: PVC (Polyvinylchloride) ●screwcap: ABS (Acrylonitrile Butadiene Styrene)

Article code	Article name	Material
Packaging		
	Tubular foil for packaging	● Film 500m x 350 mm x 40 µm ● Mixture of different Polyolefines PE (Polyethylene), LDPE (Polyethylene, low density), transparent
	Compound foil 70 µm	● PA (Polyamide) / adhesive layer / PE (Polyethylene)
	PE-Band 90 x 30 x 0.05 mm	● HDPE (high density polyethylene)
	White paper for IFU printing	● 100% wood pulp
	Roll Label 105 x 148 mm	● Paper (blank white) and Adhesive
	Adhesive tape	● BOPP (biaxially oriented polypropylene) film and rubber adhesive
	Carton 400 x 265 x 205 mm	● Cardboard
	Stretch film	● LLDPE (linear low density polyethylene)
	Stretch film foil plate	● LDPE (low density polyethylene)
	Paper MS 100 x 225 mm	● White kraft pulp 60/70/80 g/qm
Sterilization		
	Sterilization gas, Nanchang	● Ethylene oxide
	Indicator tape for EO sterilization	● Process indicator tape for Ethylene Oxide (EO) sterilization